THIS PACKAGE INSERT IS APPLICABLE FOR VETERINARY USE OUTSIDE THE US

REF	Product name	Kit configuration
983832	VET Alkaline Phosphatase (IFCC) Plus	981832 Reagent A: 15 x 4 ml 981832 Reagent B: 1 x 16 ml
983769	VET ALT / GPT (IFCC)	981769 Reagent A: 8 x 16 ml 981769 Reagent B: 8 x 4 ml
983950	VET Amylase (IFCC)	981950 Reagent A: 4 x 20 ml 981950 Reagent B: 4 x 5 ml
983771	VET AST / GOT (IFCC)	981771 Reagent A: 8 x 16 ml 981771 Reagent B: 8 x 4 ml
983888	VET Bicarbonate	981888 Reagent: 4 x 4 ml
983909	VET Bilirubin Direct	981909 Reagent A: 4 x 20 ml 981909 Reagent B: 4 x 5 ml
983952	VET Bilirubin Total (DCA)	981952 Reagent A: 4 x 20 ml 981952 Reagent B: 4 x 5 ml
983953	VET Calcium	981953 Reagent : 12 x 20 ml
983957	VET Cholesterol	981957 Reagent : 12 x 20 ml
983910	VET Cholinesterase	981910 Reagent A: 3 x 20 ml 981910 Reagent B: 3 x 5 ml
983829	VET CK (IFCC)	981829 Reagent A: 2 x 20 ml 981829 Reagent B: 2 x 5 ml
983896	VET Creatinine (Enzymatic)	981896 Reagent A: 4 x 20 ml 981896 Reagent B: 4 x 10 ml
983811	VET Creatinine (Jaffe)	981811 Reagent A: 8 x 10 ml 981811 Reagent B: 8 x 10 ml + 10 ml extra
983778	VET Gamma-GT (IFCC)	981778 Reagent A: 10 x 16 ml 981778 Reagent B: 10 x 4 ml
983780	VET Glucose (GOD-POD)	981780 Reagent : 12 x 20 ml
983779	VET Glucose (HK)	981779 Reagent A: 8 x 16 ml 981779 Reagent B: 8 x 4 ml
983955	VET HDL-Cholesterol	981955 Reagent A: 8 x 18 ml 981955 Reagent B: 8 x 6 ml
983236	VET Iron	981236 Reagent A: 10 x 20 ml 981236 Reagent B: 10 capsules 981236 Reagent C: 1 x 10 ml
983906	VET LDH (IFCC)	981906 Reagent A: 4 x 20 ml 981906 Reagent B: 4 x 5 ml
983956	VET LDL-Cholesterol	981956 Reagent A: 4 x 18 ml 981956 Reagent B: 4 x 6 ml
983846	VET Lipase (Colorimetric)	981846 Reagent A: 4 x 20 ml 981846 Reagent B: 2 x 10 ml
983905	VET Magnesium	981905 Reagent : 8 x 6 ml
983891	VET Phosphorus	981891 Reagent A: 8 x 4 ml 981891 Reagent B: 8 x 6 ml
983839	VET Pyridoxal Phosphate	981839 Reagent: 4 x 5 ml
983827	VET Total Protein Plus	981827 Reagent : 12 x 20 ml
983786	VET Triglycerides	981786 Reagent : 12 x 20 ml
983843	VET U/CSF Protein	981843 Reagent: 8 x 20 ml
983954	VET Urea	981954 Reagent A: 10 x 20 ml 981954 Reagent B: 10 x 5 ml
983788	VET Uric acid (AOX)	981788 Reagent A: 8 x 16 ml 981788 Reagent B: 8 x 4 ml

NOTE
For one assay, use always reagent components which have the same REF code and lot number.

DESCRIPTION OF THE PRODUCTS
These products are used for quantitative determinations of analytes with the Thermo Scientific™ Indiko™ clinical chemisty analyzers.

PRINCIPLE OF THE PROCEDURE
Indiko analyzer employs colorimetric end-point and kinetic, as well as turbidimetric and bichromatic reactions, with or without sample blanking.

REAGENT INFORMATION

Precautions and Preparations
Exercise the normal precautions required for handling all laboratory reagents.
The product has to be disposed of as laboratory chemical in accordance with local regulations.

Check that there are no bubbles in the bottleneck or on the surface of the reagent when you insert the reagent vials in the analyzer.
Please check the required preparation steps for the assay you are using from the table below.

Product name	Precautions	Preparation	Note
VET Alkaline Phosphatase (IFCC) Plus	The buffer is alkaline and the reagents contain sodium azide as preservative. Do not swallow. Avoid contact with skin and mucous membranes.	Add 1 ml of Reagent B into each bottle of Reagent A. Cap the vial and mix.	The p-Nitrophenol formed in the reaction is toxic. Avoid any contamination.
VET ALT / GPT (IFCC)	Contains sodium azide as preservative. Do not swallow. Avoid contact with skin and mucous membranes.	**1-reagent method (without pyridoxal-5-phosphate (P-5´-P))** Add the contents of Reagent B to the bottle of Reagent A. Cap vial and mix. **3-reagent method (with P-5´-P)** Ready for use	Simultaneous measurement of 1-reagent and 3-reagent method for the same analyte is not possible. Please choose either method and insert on-board only the reagents needed for the selected method.
VET Amylase (IFCC)	Contains sodium azide as preservative. Do not swallow. Avoid contact with skin and mucous membranes.	Ready for use	
VET AST / GOT (IFCC)	Contains sodium azide as preservative. Do not swallow. Avoid contact with skin and mucous membranes.	**1-reagent method (without pyridoxal-5-phosphate (P-5´-P))** Add the contents of Reagent B to the bottle of Reagent A. Cap vial and mix. **3-reagent method (with P-5´-P)** Ready for use	Simultaneous measurement of 1-reagent and 3-reagent method for the same analyte is not possible. Please choose either method and insert on-board only the reagents needed for the selected method.
VET Bicarbonate	Contains sodium azide as preservative. Do not swallow. Avoid contact with skin and mucous membranes.	Ready for use	Reagent is a concentrate. Water requirements for the analyzer: Use deionized or distilled water which fullfills the requirements for Clincal Laboratory Reagent Water (CLRW) of CLSI document C3-A4. The reagent is sensitive to CO_2 present in ambient air.
VET Bilirubin Direct	Reagents A and B are hazardous.	Ready for use	See separate sheet inside the kit for Hazardous- and Precautions-phrases: H290, P234, P390.
VET Bilirubin Total (DCA)	Reagents A and B are hazardous.	Ready for use	See separate sheet inside the kit for Hazardous- and Precautions-phrases: H290, H319, P280, P305+P351+P338, P337+P313
VET Calcium	Contains sodium azide as preservative. Do not swallow. Avoid contact with skin and mucous membranes. Arsenazo III is a toxic compound. Avoid contact with skin, eyes or mucous membranes.	Ready for use	Water requirements for the analyzer: Use deionized or distilled water which fulfills the requirements for Clincal Laboratory Reagent Water (CLRW) of CLSI document C3-A4.
VET Cholesterol	Contains sodium azide as preservative. Do not swallow. Avoid contact with skin and mucous membranes.	Ready for use	
VET Cholinesterase	N/A	Ready for use	
VET CK (IFCC)	Reagent A and B are hazardous. Reagents contain Imidazole < 1.5 %. The reagents contain sodium azide as preservative. Do not swallow. Avoid contact with skin and mucous membranes.	Ready for use	See separate sheets inside the kit for Hazardous- and Precautions-phrases: A: H360D, P202, P280, P308+P313. B: H315, H319, H360D, P202, P280, P302+P352, P305+P351+P338, P308+P313.
VET Creatinine (Enzymatic)	N/A	Ready for use	
VET Creatinine (Jaffe)	Reagent B is hazardous. Reagent A contains picric acid and will stain the skin.	Ready for use	See separate sheet inside the kit for Hazardous- and Precautions-phrases: H290, H314, P280, P303+P361+P353, P305+P351+P338.

thermoscientific

Product name	Precautions	Preparation	Note
VET Gamma-GT (IFCC)	The reagent contains sodium azide as preservative. Do not swallow. Avoid contact with the skin and mucous membranes.	Add 1 ml of Reagent B to each bottle of Reagent A. Cap the vial and mix.	
VET Glucose (GOD-POD)	The reagent contains sodium azide as preservative. Do not swallow. Avoid contact with skin and mucous membranes.	Ready for use	
VET Glucose (HK)	The reagents contain sodium azide as preservative. Do not swallow. Avoid contact with skin and mucous membranes.	Ready for use	
VET HDL-Cholesterol	This product contains materials of animal origin (bovine), which have been sourced from veterinary-certified healthy animals.	Ready for use	The pink instrinsic color of the reagent does not interfere with the test.
VET Iron	Reagent A is hazardous.	Dissolve *only* the contents of one capsule of Reagent B into Reagent A. Dilute Reagent C 1 + 10 with distilled water. Use the empty 20 ml bottle in the kit. Use iron free water and equipment. Check that in reagent B, the powder inside the capsule is dry and comes out easily. If not, please do not use the reagent.	See separate sheet inside the kit for Hazardous- and Precautions-phrases: H302, H315, H319, P280, P101+P312, P302+P352, P305+P351+P338.
VET LDH (IFCC)	N/A	Ready for use	
VET LDL-Cholesterol	Reagent A is hazardous. This product contains materials of animal origin (bovine), which have been sourced from veterinary-certified healthy animals.	Ready for use	See separate sheet inside the kit for Hazardous- and Precautions-phrases: Warning. H317, P261, P280, P333+P313, P362+P364
VET Lipase (Colorimetric)	Reagent B is hazardous. The reagent A contains sodium azide (NaN_3) as preservative. Do not swallow. Avoid contact with skin and mucous membranes.	Ready for use Do not mix the Lipase B reagent immediately before use. It is recommended to keep the reagent on-board over night or at the very least for a couple of hours before running a calibration and/or samples in order not to mix precipitation, if seen, with the rest of the reagent.	See separate sheet inside the kit for Hazardous- and Precautions-phrases: H319, P280, P305+P351+P338, P337+P313. As many other clinical reagents contain lipase or high concentrations of detergents, avoid carry over! Special care should be taken in combination with Triglycerides and Cholesterol reagents.
VET Magnesium	Reagent is hazardous.	Ready for use	See separate sheet inside the kit for Hazardous- and Precautions-phrases: H315, H318, P280, P302+P352, P305+P351+P338.
VET Phosphorus	Reagents A and B are hazardous.	Ready for use	See separate sheet inside the kit for Hazardous- and Precautions-phrases: H290, H314, H412, P280, P303+P361+P353, P305+P351+P338, P273.
VET Total Protein Plus	Reagent is hazardous.	Ready for use	See separate sheet inside the kit for Hazardous- and Precautions-phrases: H290, H314, P280, P303+P361+P353, P305+P351+P338.
VET Triglycerides	Contains sodium azide as preservative. Do not swallow. Avoid contact with skin and mucous membranes.	Ready for use	
VET U/CSF Protein	Reagent is hazardous.	Ready for use	See separate sheet inside the kit for Hazardous- and Precautions-phrases: H371, P260, P280, P309+P311. Reagent vials may contain precipitation. Do not mix the reagent immediately before use. It is recommended to keep the reagent on board over night or at least couple of hours before running calibration or samples.
VET Urea	Contains sodium azide as preservative. Do not swallow. Avoid contact with the skin and mucous membranes. This product contains materials of animal origin (bovine), which have been sourced from veterinary-certified healthy animals.	Ready for use	
VET Uric acid (AOX)	Contains sodium azide as preservative. Do not swallow. Avoid contact with skin and mucous membranes.	Ready for use	

Storage and Stability
The regents are stable as shown below.

Product name	Reagents, in unopened vials	Reagents, open on-board stability (days)	Note
VET Alkaline Phosphatase (IFCC) Plus	Until the expiration date printed on the label at 2…8°C	2	Keep away from light sources.
VET ALT / GPT (IFCC)	Until the expiration date printed on the label at 2…8°C	14	Keep away from light sources.
VET Amylase (IFCC)	Until the expiration date printed on the label at 2…8°C	30	
VET AST / GOT (IFCC)	Until the expiration date printed on the label at 2…8°C	14	Keep away from light sources.
VET Bicarbonate	Until the expiration date printed on the label at 2…8°C	21	
VET Bilirubin Direct	Until the expiration date printed on the label at 2…8°C	30	Keep away from light sources.
VET Bilirubin Total (DCA)	Until the expiration date printed on the label at 2…8°C	30	Do not freeze the reagents. Keep away from light sources.
VET Calcium	Until the expiration date printed on the label at 2…8°C	30	
VET Cholesterol	Until the expiration date printed on the label at 2…8°C when contamination is avoided.	30	Keep away from light sources.
VET Cholinesterase	Until the expiration date printed on the label at 2…8°C	30	Keep away from light sources.
VET CK (IFCC)	Until the expiration date printed on the label at 2…8°C	21	Keep away from light sources. If the reagents show turbidity, do not use them.
VET Creatinine (Enzymatic)	Until the expiration date printed on the label at 2…8°C	30	
VET Creatinine (Jaffe)	Until the expiration date printed on the label at 2…25°C	Reagent A: 21 Reagent B: 1	
VET Gamma-GT (IFCC)	Until the expiration date printed on the label at 2…8°C	14	Keep away from light sources.
VET Glucose (GOD-POD)	Until the expiration date printed on the label at 2…8°C	30	Keep away from light sources. De not freeze! Do not use any vial that has visible evidence of microbial growth.
VET Glucose (HK)	Until the expiration date printed on the label at 2…8°C	30	Keep away from light sources. De not freeze! Do not use any vial that has visible evidence of microbial growth.
VET HDL-Cholesterol	Until the expiration date printed on the label at 2…8°C	30	
VET Iron	Until the expiration date printed on the label at 2…8°C	21	Keep away from sunlight. Seal the plastic bag of reagent B tightly to prevent capsules to absorb moisture.
VET LDH (IFCC)	Until the expiration date printed on the label at 2…8°C	30	Keep away from light sources.
VET LDL-Cholesterol	Until the expiration date printed on the label at 2…8°C	30	
VET Lipase (Colorimetric)	Until the expiration date printed on the label at 2…8°C	30	Keep away from light sources. Do not freeze!
VET Magnesium	Until the expiration date printed on the label at 2…8°C	3	Do not use turbid reagent. Avoid contamination of the reagents. Magnesium is present in a number of products used as detergents, in tap water, etc.
VET Phosphorus	Until the expiration date printed on the label at 15…25°C	14	Reagents, once opened may be used 30 days when stored tightly sealed and protected from contamination while not in use With time, the reagent B may develop a blue color. This does not interfere with the assay. If the bias of the calibration curve is greater than 0.400 do not use the reagent.
VET Total Protein Plus	Until the expiration date printed on the label at 2…25°C	21	Keep away from light sources.

thermoscientific

Product name	Reagents, in unopened vials	Reagents, open on-board stability (days)	Note
VET Triglycerides	Until the expiration date printed on the label at 2…8°C	30	Keep away from light sources.
VET U/CSF Protein	Stable at 15…25 °C until the expiration date printed on the label when protected from light.	30	Do not use any vial that has visible evidence of microbial growth.
VET Urea	Until the expiration date printed on the label at 2…8°C	30	Keep away from light sources. Do not freeze!
VET Uric acid (AOX)	Until the expiration date printed on the label at 2…8°C	30	Keep away from light sources.

SPECIMEN COLLECTION
When processing samples in sample collection tubes, follow the instructions of the tube manufacturer carefully to avoid erroneous results. Pay special attention to the preanalytical variables such as mixing, standing time before centrifugation and centrifuge settings.

Sample Type Recommendations and Storage
Follow the general laboratory instructions when storing the samples. Biological samples should be handled and disposed of as if they were potentially infectious.

Product name	Sample type recommendations
VET Alkaline Phosphatase (IFCC) Plus	Serum or heparin plasma. Complexing anticoagulants, such as citrate, oxalate, and EDTA, must be avoided
VET ALT / GPT (IFCC)	Clear serum or heparin plasma
VET Amylase (IFCC)	Serum, heparin plasma or urine If the urine sample is tested with the delay its pH should be adjusted to approximately 7.0 Only freeze once! Saliva contaminates the sample. Discard contaminated specimens!
VET AST / GOT (IFCC)	Unhemolyzed serum or heparin plasma
VET Bicarbonate	Serum or plasma (Li-heparin) Keep sample tubes closed.
VET Bilirubin Direct	Clear, non hemolyzed serum or Li-heparin plasma Bilirubin is unstable in the samples. Protect the samples from direct solar or artificial light.
VET Bilirubin Total (DCA)	Serum or Li-heparin plasma Bilirubin is unstable in the samples. Protect the samples from direct solar or artificial light.
VET Calcium	Serum, Li-heparin plasma or urine Do not use citrate, oxalate or EDTA because they remove calcium by complex formation. Any substance which either chelates or contains calcium will interfere with the test. Urine samples should be acidified to prevent calcium salt precipitation. 24-hour specimens may be collected in 20 to 30 mL of 6 N HCl. Random specimens may be acidified with 1 - 2 mL of 6 M HCl
VET Cholesterol	Serum or heparin plasma
VET Cholinesterase	Serum or heparin plasma
VET CK (IFCC)	Unhemolyzed serum or heparin plasma Hemolysis can lead to deviating results. Avoid exposure of samples to strong light.
VET Creatinine (Enzymatic)	Serum, Li-heparin plasma and urine
VET Creatinine (Jaffe)	Serum, Li-heparin plasma and urine Avoid using hemolyzed specimen. Urine samples must be diluted with distilled or deionized water according to dilution schema desired (e.g. 1+50), done automatically on Indiko analyzer.
VET Gamma-GT (IFCC)	Serum or heparin plasma Fluoride, oxalate or citrate should not be used because of the inhibition of Gamma-GT activity
VET Glucose (GOD-POD)	Unhemolyzed serum or heparin plasma The serum and plasma samples should be separated from the cells as soon as possible after collection in order to avoid glycolysis. If the serum or plasma sample is not separated or analyzed without delay a glycolytic inhibitor should be used.
VET Glucose (HK)	Unhemolyzed serum, heparin plasma or urine The serum and plasma samples should be separated from the cells as soon as possible after collection in order to avoid glycolysis. If the serum or plasma sample is not separated or analyzed without delay a glycolytic inhibitor should be used.
VET HDL-Cholesterol	Serum or heparin plasma Specimens should preferably be analyzed on the day of collection. Samples containing precipitate should be centrifuged before performing the assay
VET Iron	Serum or heparin plasma Avoid using hemolyzed and lipemic specimens. Do not use EDTA, citrate or oxalate plasma.
VET LDH (IFCC)	Non-hemolyzed serum or heparin plasma Be careful with samples as they may be contaminated with platelets, which contain high concentrations of LDH. The samples should be centrifuged properly.
VET LDL-Cholesterol	Fasting and non-fasting serum and heparin plasma
VET Lipase (Colorimetric)	Serum or heparin plasma
VET Magnesium	Non hemolyzed serum or heparin plasma separated as soon as possible from the cells Do not use metal complexing anticoagulants like EDTA
VET Phosphorus	Serum, Li-heparin plasma or urine Do not use anticoagulants such as citrate, oxalate or EDTA because they interfere with the formation of the phosphomolybdate complex. Use clear, unhemolyzed serum or plasma, separated from the erythrocytes as soon as possible. Erythrocytes contain organic phosphates which can hydrolyze on standing or can be enzymatically cleaved by phosphatases. Inorganic phosphates can then leak through the cell walls, increasing the concentration in serum or plasma. Urine samples should be fresh and acidified to avoid precipitation of phosphate complexes. 24-hour specimens may be collected in 20 to 30 ml of 6 N HCl. Random specimens may be acidified with 1-2 ml of 6 N HCl (3).
VET Total Protein Plus	Serum or heparin plasma Unhemolyzed sample is preferred for analysis. Plasma gives approximately 1 - 5 g/l higher values due to fibrinogen present.
VET Triglycerides	Serum or heparin plasma The specimen should be collected after a 10 - 14 hour fast.
VET U/CSF Protein	Urine or cerebrospinal fluid Blood contamination should be avoided when obtaining the cerebrospinal fluid sample, because it may give false positive results. The sample must be centrifuged immediately.
VET Urea	Serum, heparin plasma or urine Do NOT use serum preserved with fluoride. Ammonium ions interfere; therefore, do not use ammonium heparin as anticoagulant for collection of plasma. Urine samples must be diluted with distilled or deionized water according to dilution schema desired (e.g. 1+50), done automatically on analyzers.
VET Uric acid (AOX)	Serum, heparin plasma or urine Urine samples must be diluted with distilled or deionized water according to dilution schema desired (e.g. 1+10), done automatically on analyzers. Add NaOH to make urine alkaline.

TEST PROCEDURE
To enable the VET testing, either create user-definable tests or modify applications available in the analyzer system database. Applications in the analyzer database must not be used as such. Creating user-definable tests or modifying existing tests requires using a dongle that enables saving these tests. The dongle is installed in the Indiko analyzer's computer's USB port.

To install the dongle in the computer's USB port:
1. If the computer has a used dongle, remove it from the USB-port.
2. Install a new dongle to the USB-port in the keyboard.
3. In the Indiko analyzer user interface select **F5 > Configuration > Additional**.
4. Check the number of user-definable tests from the display.

NOTE: The dongle is computer-specific. When it is once connected to certain computer it cannot be connected to another one.

To create a test procedure:
In the Indiko analyzer user interface
1. Select **F3 > Reagent definition** to add reagents.
2. Select **F4 > Cal/QC selection > QC profile** definition to add QC profiles.
3. Select **F4 > Cal/Ctrl definition** to add controls and calibrators.
4. Select **F3 > Test definition** to define the test parameters.
5. Select **F4 > Cal/Ctrl definition > Lot concentrations** to define the lot concentrations for controls and calibrators.

thermoscientific

To modify the existing test in the database:
In the Indiko analyzer user interface

1. Select **F3 > Test definition**.
2. Select the system test to be edited and click **New**.
3. Select **Yes** when the system asks 'Do you want to copy parameters of current test.'
4. Give a name to the test and click **Confirm**.
5. Edit the test parameters as needed and click **Save** to save the changes in each view.
6. Take the new test into use by changing the In use setting to **Yes** in the Test definition > **Info** view. Click **Save**.

NOTE: Each laboratory is responsible for the application validation with the used sample matrix to ensure the test performance.

Materials available but not provided
Please ensure that the below mentioned calibrators and controls are suitable for the modified applications. Enter the values found in the value sheets manually. Calibrator and control values must be verified for the used applications.

Product name	Calibrator	Quality controls	Other products
VET Alkaline Phosphatase (IFCC) Plus	VET eCal	VET Nortrol, VET Abtrol	
VET ALT / GPT (IFCC)	VET eCal	VET Nortrol, VET Abtrol	Pyridoxal-5'-phosphate for 3-reagent method
VET Amylase (IFCC)	VET eCal	Serum/ plasma: VET Nortrol, VET Abtrol Urine: VET uTrol, VET uTrol High	
VET AST / GOT (IFCC)	VET eCal	VET Nortrol, VET Abtrol	Pyridoxal-5'-phosphate for 3-reagent method
VET Bicarbonate	VET sCal	VET Nortrol, VET Abtrol	
VET Bilirubin Direct	VET sCal	VET Nortrol, VET Abtrol	
VET Bilirubin Total (DCA)	VET sCal	VET Nortrol, VET Abtrol	
VET Calcium	VET sCal	Serum/ plasma: VET Nortrol, VET Abtrol Urine: VET uTrol, VET uTrol High	
VET Cholesterol	VET sCal	VET Lipotrol, VET Nortrol, VET Abtrol	
VET Cholinesterase	VET sCal	VET Nortrol, VET Abtrol	
VET CK (IFCC)	VET eCal	VET Nortrol, VET Abtrol	
VET Creatinine (Enzymatic)	VET sCal	Serum/ plasma: VET Nortrol, VET Abtrol Urine: VET uTrol, VET uTrol High	Washing liquid, WashFluid, to eliminate risks of cross-contamination
VET Creatinine (Jaffe)	VET sCal	Serum/ plasma: VET Nortrol, VET Abtrol Urine: VET uTrol, VET uTrol High	
VET Gamma-GT (IFCC)	VET eCal	VET Nortrol, VET Abtrol	
VET Glucose (GOD-POD)	VET sCal	VET Nortrol, VET Abtrol	
VET Glucose (HK)	VET sCal	Serum/ plasma: VET Nortrol, VET Abtrol Urine: VET uTrol, VET uTrol High	
VET HDL-Cholesterol	VET HDL/LDL Calibrator	VET Lipotrol, VET Lipotrol HDL/LDL Abnormal	
VET Iron	VET sCal	VET Nortrol, VET Abtrol	
VET LDH (IFCC)	VET eCal	VET Nortrol, VET Abtrol	
VET LDL-Cholesterol	VET HDL/LDL Calibrator	VET Lipotrol, VET Lipotrol HDL/LDL Abnormal	Washing liquid, WashFluid, to eliminate risks of cross-contamination
VET Lipase (Colorimetric)	VET eCal	VET Nortrol, VET Abtrol	
VET Magnesium	VET sCal	VET Nortrol, VET Abtrol	
VET Phosphorus	VET sCal	Serum/ plasma: VET Nortrol, VET Abtrol Urine: VET uTrol, VET uTrol High	
VET Total Protein Plus	VET sCal	VET Nortrol, VET Abtrol	
VET Triglycerides	VET sCal	VET Lipotrol, VET Nortrol, Abtrol	
VET U/CSF Protein	VET sCal	VET uTrol High	Washing liquid, WashFluid, to eliminate risks of cross-contamination
VET Urea	VET sCal	Serum/ plasma: VET Nortrol, VET Abtrol Urine: VET uTrol, VET uTrol High	
VET Uric acid (AOX)	VET sCal	Serum: VET Nortrol, VET Abtrol Urine: VET uTrol, VET uTrol High	

Calibration
Use calibrators VET sCal (code 983831), VET eCal (code 983830), and VET HDL/LDL Calibrator (code 983657) according to the instructions provided for your analyzer.

Quality Control
Use quality control (QC) samples at least once a day and after each calibration and every time a new bottle of reagent is used. It is recommended to use at least two levels (low and high) of controls or sample pools. Always follow the local, state and federal regulations in performing QC.

Codes of available controls:
VET Nortrol, code 983043
VET Abtrol, code 983044
VET Lipotrol, code 983653
VET Lipotrol HDL/LDL Abnormal, code 983907
VET uTrol, code 983821
VET uTrol High, code 983822

The Control intervals and limits must be adapted to the individual laboratory requirements. The results of the quality control samples should fall within the limits pre-set by the laboratory.

Calculation of results
The results are calculated automatically by the analyzer according to the application using the calibration curve.

LIMITATIONS OF THE PROCEDURE
Assays to be used with companion animal samples, not with livestock samples. The manufacturer has not validated the performance, therefore each laboratory must validate the performance with the samples used.

MEASURING RANGE

Product name	Primary	Secondary
VET Alkaline Phosphatase (IFCC) Plus	20 – 1000 U/l (0.33 - 16.7 µkat/l)	20 – 10 000 U/l (0.33 - 167 µkat/l)
VET ALT / GPT (IFCC)	1-reagent method: 4 - 200 U/l (0.07 - 3.33 µkat/l) 3-reagent method: 5 - 250 U/l (0.08 - 4.17 µkat/l)	1-reagent method: 4 - 1200 U/l (0.07 - 20.00 µkat/l) 3-reagent method: 5 - 1500 U/l (0.08 - 25.00 µkat/l)
VET Amylase (IFCC)	10 - 1500 U/l (0.17 – 25 µkat/l)	10 - 7500 U/l (0.17 - 125 µkat/l)
VET AST / GOT (IFCC)	1-reagent method: 4 - 350 U/l (0.07 - 5.83 µkat/l) 3-reagent method: 4 - 300 U/l (0.07 - 5.00 µkat/l)	1-reagent method: 4 - 2100 U/l (0.07 – 35.0 µkat/l) 3-reagent method: 4 - 1800 U/l (0.07 - 30.0 µkat/l)
VET Bicarbonate	5 – 40 mmol/l (5 – 40 mEq/l)	Please dilute samples over test limit high manually 1+1 using fresh deionized or distilled water.
VET Bilirubin Direct	1.0 - 170 µmol/l (0.06 – 9.9 mg/dl)	1.0 - 1190 µmol/l (0.06 – 69.6 mg/dl)
VET Bilirubin Total (DCA)	1.0 - 400 µmol/l (0.06 - 23 mg/dl)	1.0 - 2000 µmol/l (0.06 - 117 mg/dl)
VET Calcium	Serum: 0.4 - 5.0 mmol/l (1.6- 20.1 mg/dl) Urine: 0.4 - 5.0 mmol/l (1.6- 20.1 mg/dl)	Serum: N/A Urine: 0.4 - 12 mmol/l (1.6 – 48.1 mg/dl)
VET Cholesterol	0.2 - 20.0 mmol/l (8 - 773 mg/dl)	0.2 - 45.0 mmol/l (8 - 1737 mg/dl)
VET Cholinesterase	400* - 16000* U/l (6.6* – 266.7* µkat/l)	100 – 90000* U/l (1.7 – 1500* µkat/l) *The values are related to the concentration of the highest calibrator and are lot dependent.
VET CK (IFCC)	10 – 1000 U/l (0.17 – 17 µkat/l)	10 – 9000 U/l (0.17 – 150 µkat/l)
VET Creatinine (Enzymatic)	Serum, plasma: 10 - 2500 µmol/l (0.11 - 28 mg/dl) Urine: 0.2 - 40 mmol/l (2.3 - 452 mg/dl)	Serum, plasma: 10 - 10000 µmol/l (0.11 - 113 mg/dl) Urine: 0.01 - 200 mmol/l (0.11 - 2260 mg/dl)
VET Creatinine (Jaffe)	Serum, plasma: 20 - 2500 µmol/l (0.2 - 28.3 mg/dl) Urine: 0.5 - 40.8 mmol/l (5.7 - 462 mg/dl)	Serum, plasma: 20 - 10000 µmol/l (0.2 - 113.0 mg/dl) Urine: 0.01 - 80.8 mmol/l (0.1 - 914 mg/dl)
VET Gamma-GT (IFCC)	15 - 600 U/l (0.25 - 10 µkat/l)	15 - 6000 U/l (0.25 - 100 µkat/l)
VET Glucose (GOD-POD)	0.3 - 20.0 mmol/l (5 - 360 mg/dl)	0.3 - 60.0 mmol/l (5 - 1080 mg/dl)
VET Glucose (HK)	0.3 - 40.0 mmol/l (5 - 720 mg/dl)	0.3 - 120.0 mmol/l (5 - 2160 mg/dl)
VET HDL-Cholesterol	0.10 - 3.88 mmol/l	0.10 - 11.64 mmol/l
VET Iron	1.0 - 180.0 µmol/l (6 - 1006 µg/dl)	1.0 - 540.0 µmol/l (6 - 3017 µg/dl)
VET LDH (IFCC)	10 - 600 U/l (0.17 - 10.0 µkat/l)	10 - 3000 U/l (0.17 - 50.0 µkat/l)

thermoscientific

Product name	Primary	Secondary
VET LDL-Cholesterol	0.10 - 11.0 mmol/l (4 - 425 mg/dl)	0.10 - 33.0 mmol/l (4 - 1276 mg/dl)
VET Lipase (Colorimetric)	16 - 300 U/l	16 - 1500 U/l
VET Magnesium	0.1 - 1.7 mmol/l (0.2 - 4.1 mg/dl)	0.1 – 5.1 mmol/l (0.2 – 12.4 mg/dl)
VET Phosphorus	Serum, plasma: 0.08 - 4.0 mmol/l (0.2 - 12.4 mg/dl) Urine: 2.5 - 80.0 mmol/l (7.7 - 247.8 mg/dl)	Serum, plasma: 0.08 - 20.0 mmol/l (0.2 - 61.9 mg/dl) Urine: 0.1 - 194.0 mmol/l (0.3 - 600.8 mg/dl)
VET Total Protein Plus	3 - 100 g/l (0.3 – 10.0 g/dl)	3 - 250 g/l (0.3 – 25.0 g/dl)
VET Triglycerides	0.05 - 11.00 mmol/l (4 - 973 mg/dl)	0.05 - 55.00 mmol/l (4 - 4867 mg/dl)
VET U/CSF Protein	0.080 - 3.000 g/l (8.0 - 300 mg/dl)	0.080 - 27.00 g/l (8.0 - 2700 mg/dl)
VET Urea	Serum, plasma: 1.0 - 25.0 mmol/l (6 - 150 mg/dl) Urine: 51.0 - 1275 mmol/l (306 – 7650 mg/dl)	Serum, plasma: 1.0 - 75.0 mmol/l (6 - 450 mg/dl) Urine: 4.0 - 2525 mmol/l (24 – 15150 mg/dl)
VET Uric acid (AOX)	Serum, plasma: 10 - 1200 µmol/l (0.2 - 20.0 mg/dl) Urine: 0.60 - 13.2 mmol/l (10.2 - 224 mg/dl)	Serum, plasma: 10 - 7200 µmol /l (0.2 - 120.0 mg/dl) Urine: 0.01 - 61.2 mmol /l (0.2 - 1028 mg/dl)

WASTE MANAGEMENT
Please refer to local legal requirements. It is recommended to empty the analyzer cuvette waste bin and wastewater daily. Emptying should be done immediately after the analysis when using hazardous reagents/solutions.
Note: If using reagents/solutions that react with each other, cuvette waste bin and wastewater should be emptied and washed between use of these reagents.

MANUFACTURER

Thermo Fisher Scientific Oy
Ratastie 2, FI-01620, Vantaa, Finland

edfu.thermofisher.com

CONTACT INFORMATION
Postal address: Thermo Fisher Scientific Oy,
Ratastie 2, P.O. Box 100, FI-01621 Vantaa, Finland
Tel. +358 10 329 200
www.thermofisher.com

NOTE
Report any serious incident that happens in relation to the device to the manufacturer and when required by law to the regulatory authority of the country in which the user is established.

Date of revision (yyyy-mm-dd)
2024-03-25

Changes from previous revision
The line beside the text indicates changes.

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Postadress: Thermo Fisher Scientific Oy,