
Large menu of gel electrophoresis *on HYDRASYS 2 SCAN FOCUSING*

Semi-Automated Multi-Assay Analyzer Proven Technology





sebia

The new language of life

Main *Strengths*

HYDRASYS 2 SCAN FOCUSING gel electrophoresis instrument provides your laboratory with high-quality results.

Combining the HYDRASYS 2 SCAN FOCUSING with our front-end pipetter, ASSIST, improves worflow and traceability.



High Quality Result

- Have confidence when reporting your patients' results with gel electrophoresis assays.

Multi-Functional Analyzer

- Ability to migrate, stain, & scan gels on one system.



Powerful Software Solution

- PHORESIS Core Software provides tools for interpretation, data & result management and networking.

Workflow Traceability

- Automation and standardization of samples preparation with ASSIST.


All-in-One System Gel Electrophoresis

- All-in-one system performing migration, staining & scanning
- Small footprint fits any size work area
- Intuitive touchscreen interface


Gel Automation with the ASSIST


- Preparation and automatic dispensing of samples from primary tube and antisera for dynamic mask application
- Auto-dilutes patient samples for IFE and CSF testing
- Traceability of primary tubes with positive sample ID with PHORESIS

Multi gel sizes
to fit your Laboratory's needs




Protein (Serum)
PROTEIN(E) and β1-β2
Gels sizes - 7, 15, 30, 54






Immunofixation (Serum/Urine)
Immunofixation
Gels sizes - 1, 2, 4, 9


BENCE JONES
Gels sizes - 1, 2, 4




HYDRASHIFT
daratumumab/isatuximab
Immunofixation
Gels sizes - 2, 4




URINE PROFIL(E)
Immunofixation
Gels sizes - 2, 4






Iso-Electro Focusing - CSF & A1AT
CSF ISOFOCUSING
Gels sizes - 3, 9

Alpha 1 - Antitrypsin
Gels sizes - 18




von WILLEBRAND MULTIMERS
Immunofixation
Gels sizes - 5,11




Lipids
LIPOPROTEIN(E) - LIPO + Lp(a)
Gels sizes - 7, 15, 30





Hemoglobin
HEMOGLOBIN(E) ACID(E) & Alkaline
Gels sizes - 7, 15



Enzymes
Isoforms: Lactate Dehydrogenase (LD), Creatine Kinase (CK) & Phosphatase Isoenzymes (ISO-PAL)
Gels sizes - 7, 15

LARGE RANGE OF ASSAYS

- Serum & Urine Protein Electrophoresis (SPE/UPE)
- Serum & Urine Immunofixation (IFE, BENCE JONES)
- HYDRASHIFT daratumumab/isatuximab
- CSF and Alpha-1 Antitrypsin by ISOFOCUSING
- von WILLEBRAND MULTIMERS
- Iso-Enzymes, B2 transferrin,
- LDL/HDL cholesterol, lipoproteins
- HEMOGLOBIN(E) Electrophoresis

Ordering *information*

Instruments

- HYDRASYS 2 SCAN FOCUSING 1202
- ASSIST 1218

Reagents*

- HYDRAGEL 7/15 PROTEIN(E) 4100/4120
- HYDRAGEL 30/54 PROTEIN(E) 4140/4145
- HYDRAGEL 7/15 β1-β2 ..... 4101/4121
- HYDRAGEL 30/54 β1-β2 ..... 4141/4146
- HYDRAGEL 7 HR 4102
- HYDRAGEL 15 HR 4122
- HYDRAGEL 5 PROTEINURIE 4115
- HYDRAGEL 6 IF Penta 4341
- HYDRAGEL 12 IF Penta 4342
- HYDRAGEL 1/2 IF 4801/4802
- HYDRAGEL 4/9 IF 4804/4809
- HYDRAGEL 2 IF/BJ (HR) 4803
- HYDRAGEL 1 BENCE JONES 4821
- HYDRAGEL 2 BENCE JONES 4822
- HYDRAGEL 4 BENCE JONES 4824
- HYDRAGEL 2 URINE PROFIL(E) 4831
- HYDRAGEL 4 URINE PROFIL(E) 4832
- HYDRASHIFT 2/4 daratumumab (20 TESTS) 4639
- HYDRASHIFT 2/4 isatuximab (40 TESTS) 4642

Reagents*

- HYDRAGEL 5 von WILLEBRAND MULTIMERS 4359
- HYDRAGEL 11 von WILLEBRAND MULTIMERS 4360
- HYDRAGEL 3 CSF ISOFOCUSING 4353
- HYDRAGEL 9 CSF ISOFOCUSING 4355
- HYDRAGEL 18 A1AT ISOFOCUSING 4357
- HYDRAGEL 15 LIPO + Lp(a) 4124
- HYDRAGEL 30 LIPO + Lp(a) 4144
- HYDRAGEL 15 LIPOPROTEIN(E) 4134
- HYDRAGEL 30 LIPOPROTEIN(E) 4138
- HYDRAGEL 15 LDL/HDL CHOL Direct 4135
- HYDRAGEL 30 LDL/HDL CHOL Direct 4139
- HYDRAGEL 7 ISO-LDH 4110
- HYDRAGEL 15 ISO-LDH 4130
- HYDRAGEL 30 ISO-LDH 4136
- HYDRAGEL 7 ISO-CK 4111
- HYDRAGEL 15 ISO-CK 4131
- HYDRAGEL 30 ISO-CK 4137
- HYDRAGEL 7 ISO-PAL 4112
- HYDRAGEL 15 ISO-PAL 4132
- HYDRAGEL 6 β2 TRANSFERRIN(E) 4853

Scan here **for more information on** HYDRASYS 2 SCAN FOCUSING

* Non-exhaustive list, contact your local representative for more information about the Sebia gel range.

For further information, please find your Sebia local representative to contact on our Website **www.sebia.com/sebia/worldwide-presence/.**

CE

Sebia | RCS Evry 672 041 902 | Parc technologique Léonard de Vinci | 27, rue Léonard de Vinci | CP 8010 Lisses | 91008 Evry cedex | France | Tel.: +33 1 69 89 80 80 | E-mail: sebia@sebia.com | **www.sebia.com**

©2025 **Sebia** - RCS Evry 672 041 902 | Parc technologique Léonard de Vinci | 27, rue Léonard de Vinci | CP 8010 Lisses | 91008 Evry cedex | France - Ref.: H2SB1EN - 04/2023
Non-contractual illustration - Sebia reserves the right to modify at any time and without notice the information contained in this document intended for healthcare professionals. This document contains information intended for wide distribution and may therefore contain product details or information that is not available or valid in your country. Carefully read the instructions in the reagent package inserts and instrument manuals.